

DDDDK-Tag(binds to flag sequnence) rabbit pAb

Catalog No: YG0004

Reactivity: Species independent

Applications: WB;ELISA;IP;IF

Target: FLAG-tag

Gene Name: Flag tag; Flag-tag, DDDDK TAG, DDDDK-TAG, DYKDDDDK tag, DYKDDDDK-

tag

Immunogen: DDDDK synthetic peptide conjugated to KLH.

Specificity: FLAG-tag Polyclonal Antibody detects FLAG-tagged recombinant proteins or

FLAG-tagged proteins overexpressed in cells.

Formulation : Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Source: Polyclonal, Rabbit, IgG

Dilution: WB 1:1000 - 1:3000. ELISA: 1:5000-20000.. IF 1:50-200

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Concentration: 1 mg/ml

Storage Stability: -15°C to -25°C/1 year(Do not lower than -25°C)

Background: The DYKDDDDK (FLAG) peptide has been used extensively as a general tag in

expression vectors. This peptide can be expressed and detected with the protein of interest as an amino-terminal or carboxy-terminal fusion. N-terminal FLAG vectors provide an Ek cleavage site for removal of the fusion tag. The FLAG peptide is likely to be located on the surface of a fusion protein because of its hydrophilic nature. As a result, the FLAG peptide is more likely to be accessible to antibodies. A FLAG-tag can be used in many different assays that require recognition by an antibody, such as western blotting, immunocytochemistry,

immunoprecipitation, flow cytometry, protein purification, and in the study of protein-protein interactions, cell ultrastructure, and protein localization and so on.

1/2



Tag: orthogonal,overexpres,ip


Sort : 603

No4: 1

Host: Rabbit

Modifications: Unmodified

Products Images



Western Blot analysis using FLAG-tag Polyclonal Antibody against HEK293 cells transfected with vector overexpressing FLAG tag (1) and untransfected (2). Secondary antibody(catalog#:RS0002) was diluted at 1:20000